Long read sequencing: High-fidelity amplification of disease-associated complex genomic regions using DeCodi-Fipolymerase.

J. Hsieh, J. Gimpel, J. Sáez, J. Cáceres, C. Valdivia, I. Arriagada, T. Akentjew, S. Kamkar, I. Marx, M. Rozas Kura Biotech Spa, Puerto Varas, Región de Los Lagos, Chile

Background

Repeat expansion disorders such as Amyotrophic Lateral Sclerosis (ALS), Friedreich's Ataxia (FrA), Fragile X (FrX) and Huntington's Disease (HD) are notoriously difficult to

Repeat Expansion Samples: DNA from unaffected control and affected individuals (Coriell Institute) were used to evaluate amplification performance across **expanded repeat loci** for each disease by different polymerases. **Table 1** summarizes target

Methods

Table 1: Summary of genes, repeat motifs, template DNA, repeat counts, and expected amplicon sizes for each disease studied.

Poster No. #B-224

KURA

detect with short-read sequencing (fragments shorter than 300 bp)	genes, repeat motifs, template DNA sources, repeat counts and expected product sizes for each disease.	Disease	Gene	Repeat Motif	Template DNA (Coriell Sample)	Repeat Count (Allele 1/Allele 2)	Expected Size (Allele 1/Allele 2)
 Amplifying GC-rich, long tandem repeats often results in polymerase slippage or dropout, which can compromise 	Long-Range PCR Samples: To evaluate amplification performance on large genomic	Amvotrophic			Unaffected control	< 20	143 bp
diagnostic accuracy when an amplification step is	DeCodi-Fi's ability to amplify long human DNA fragments compared to other	BRCAI (8,4 kb, 23,6 kb and 36,4 kb) and SMNI (28.2 kb), assessing to amplify long human DNA fragments compared to other Lateral		GGGGCC	ALS (ND12780)	22/70	178/466 bp
 Long-read sequencing platforms (e.g., Nanopore) 	high-fidelity polymerases.	Scierosis (ALS)			ALS (ND12667)	5/~950	76/5746 bp
provide a more direct and reliable approach for resolving repeat expansions, but their success depends	Polymerases tested: DeCodi-Fi (Kura Biotech), Platinum SuperFi II (Thermo Fisher Scientific) and Q5 (NEB). Relevant reaction buffers for each enzyme were used	Friedreich's Ataxia (FrA)	X25	GAA	Unaffected control	< 33	520 bp
heavily on the performance of the high-fidelity	according to manufacturer's instructions. Where additional enhancers to PCR reactions Ataxia were used, the same additions were made to competitor reactions.				FrA (NA16243)	670/1117	2510/4010 bp
templates such as long, GC-rich regions (e.g. 6.5 kb with	Cycling conditions: Cycling protocols were customized and optimized for specific	nditions: Cycling protocols were customized and optimized for specific s and targets, using techniques such as touchdown PCR, controlled ramp wobble" temperatures. The same amplification protocol was applied across		Unaffected control	<44	350 bp	
	primer pairs and targets, using techniques such as touchdown PCR, controlled ramp rates and "wobble" temperatures. The same amplification protocol was applied across			000	FrX (NA7862)	> 200-230	2339bp/>3239bp
Objective	all polymerases for each target.				Unaffected control	< 36	105 bp
To evaluate the performance of the high-fidelity polymerase DeCodi-Fi compared to other high-fidelity	Detection method: Results were visualized by agarose gel electrophoresis. Long-read	Disease (HD) <i>IT15</i> CAC	CAG	HD (NA20253)	22/100	119/353 bp	
enzymes in amplifying complex, long, and extremely GC-rich repeat regions.	sequencing was performed on the Fromethion (Oxford Nahopore rechnologies).				HD (NA13509)	15/70	98/263 bp

Results: Amplifying the Toughest Targets (6.5 kb GC Repeats)

Figure 1. Amplification of disease-associated repeat expansions using three high-fidelity polymerases. Agarose gel electrophoresis showing PCR amplification of repeat-containing loci from unaffected control and disease-affected samples across four disorders: ALS (A), FrA (B), FrX (C), and HD (D). Three high-fidelity polymerases; DeCodi-Fi (DC), Platinum SuperFi II (S), and Q5 (Q) were tested under equivalent conditions. For unaffected control amplification, DeCodi-Fi and Q5 were used with their respective GC buffer/enhancer, and SuperFi with its reaction buffer. Additional enhancers were included to amplify disease-expanded alleles. Light blue arrows indicate detected PCR products for each allele. Some gels show primer dimers or non-specific bands due to pre-cleanup loading, intended to enable raw product comparison across conditions. M: 100 bp ladder; L: 1 kb Ladder.

A. Amyotrophic	c Lateral S	clerosis (ALS)		B. Friedreich's A	Ataxia (FrA)	C. Fragile X Syndro	me (FrX)		D. Huntington's Disease (HD)						
DeCodi-Fi	Control	ND12780 ND12	667	DeCodi-Fi amplified	Control NA16243	DeCodi-Fi was the only	Control FrX	NA7862			Control		HD NA20253	HD NA135	509
uniquely amplified	L DC S Q	DC S Q M DC S	QL	both alleles in the	M DC S Q L DC S Q	enzyme to amplify the	M DC S Q L	DC S Q	amplified both	М	DC S	Q	M DC S	Q M DC	S Q
both alleles in				FrA sample (2510		expanded 6500 bp	*		affected and control				- 100 mil		and with
ND12780 (178 bp /		¢	5746 bp	bp / 4010 bp).		mosaic allele (100%			samples (98–263						
466 bp) and	-			SuperFi and Q5		GC), which is often	=	6500 bp	bn) DeCodi-Fi			- * T			
ND12667 (75 bp /		The second second second		produced weaker	E 4010 bp	missed and present		<	vielded the cleanest						
5746 bp). SuperFi	1000	COLUMN AND ADDRESS	-	bands. All enzymes	2510 bp	only in a	• •		hands Products						
and Q5 failed to	-	-	-	successfully		subpopulation of cells	_	< 2339 bp	were sequenced via						
amplify specific		1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1	100 C 100 C	amplified the	٤,	in this FrX affected			Oxford Nanonore						
bands, producina		-	-	control (~520 bp).	-	sample. The 2339 bp	-			_					
only smears or			100		-	allele and the	-					1.00	The second second		



MAIN FINDING:

DeCodi-Fi robustly amplified 100% GC targets up to ~500 bp and was the only enzyme to detect larger repeats (2000–6500 bp), outperforming competitors.

Oxford Nanopore Sequencing

Amplifying the Longest Targets

Huntington's disease alleles were selected for comparing sequencing results between polymerases since it was the only disease where both unaffected and affected samples could be amplified by all polymerases. DeCodi-Fi generally achieved higher coverage and in some cases more spanning reads than Platinum SuperFi and Q5.

Table 2: Motif counts observed and sequencing metrics for unaffected individual and Huntington's disease (HD) patients for Allele1(A1)/Allele2(A2). DeCodi-Fi (DC), Platinum SuperFi II (S), and Q5 (Q) polymerases were compared. The median motif count refers to the number of motifs detected in the sequencing results for each allele. Coverage corresponds to the total number of sequencing reads aligned over the repeat region, and spanning reads indicate the number of reads that fully span the repeat motif.

DeCodi-Fi successfully amplified human DNA fragments up to 36.4 kb, while other polymerases were limited to amplifying fragments no larger than 8.4 kb.

	BRC	A1 8.4	kb		BRC	AI 23.6	kb	BRCA136.4			kb	BRCA12			28.2 kb	
L	DC	S	Q	L	DC	S	Q	L	DC	S	Q	L	DC	S	Q	
		277					an a		-					100		
											2					
								-	_		Sec.					
					-								-			

Conclusion

• DeCodi-Fi reliably amplified all four disease-associated repeat expansions, including long and highly GC-rich alleles.

- It was the only polymerase to amplify the 6.5 kb, 100% GC ([CGG]>1000) repeat in Fragile X and the 5.7 kb GGGGCC expansion in ALS.
- In long-range PCR, DeCodi-Fi amplified human genomic targets up to 36.4 kb, whereas other polymerases were limited to 8.4 kb.

Genome	Expected motif count (A1/A2)	Enzyme	Median motif count(A1/A2)	Coverage	Spanning reads (A1/A2)
		DC	17/17	9958	2367/2067
Unaffected	<36	S	18/18	9948	2104/2254
		Q	18/18	8122	1856/1675
	22/100	DC	22/103	8930	2652/230
HD NA20253		S	22/106	8016	2526/380
		Q	22/104	5706	1766/279
		DC	15/71	8568	3404/511
HD NA13509	15/70	S	15/71	748	252/93
		Q	22/71	4720	2527/381



Figure 2. Long-range amplification of clinically relevant targets, including Breast Cancer I gene (BRCAI) and Survival of Motor Neuron I gene (SMNI). DeCodi-Fi (DC) successfully amplified all four large fragments (BRCA1 8.4 kb, 23.6 kb, 36.4 kb, and SMN1 28.2 kb), producing strong, specific bands. This capability is critical for long-read sequencing, where full-length amplification enables comprehensive analysis of structural variants and mutations. Platinum SuperFi II (S) amplified only the 8.4 kb fragment, while Q5 (Q) produced a faint band at 8.4 kb and both failed to amplify larger targets. L: 1 kb ladder.

• Sequencing of Huntington's disease samples showed higher coverage and more spanning reads with DeCodi-Fi.

• These results support DeCodi-Fi as a robust tool for long-read sequencing of complex genomic regions.